

VN415

Feature Overview



Interacoustics[®]

leading diagnostic solutions

Connect and Go Testing

Start the software, place the goggles on the patient and you are ready to test. The ease of use of the VN415 is a major factor in its success. There is no bulky equipment to get in the way.

High Quality Video and Analysis

The VN415 performs real time analysis. Rapid data transfer (105 Hz binocular, 174 Hz monocular) allows much more data to be analyzed, giving very precise measurements and outstanding video quality.

Vestibular Test Battery

The VN415 test battery includes Spontaneous Nystagmus, Dix-Hallpike, Position and Bithermal Caloric tests. You can design your own test protocols and load these with a single mouse click. And there is a hand held infra red device that allows you to remotely step through tests.

Remote Controlled Tests

It is possible to step through the tests using a remote device. This can either be a USB connected footswitch or an infrared hand held device.

Sharing Information

VN415 has export functions designed around referral and consultation with balance testing colleagues. Results as well as recordings can be exported. Colored printed reports are available in each test.

Electronic Medical Resources / Networking (EMR)

Export to PDF or XML

Automatic and Manual Analysis

The system produces real time accurate measurements with the user integrated editing functions