

RhinoCycle

Feature Overview

Standalone unit

RhinoCycle is a standalone unit, and does not work with the SRE2000.

Easy to use software for analysis of monitored data on own pc

A complete picture of the nasal cycle is monitored, and data are imported to a pc, for further analysis. At termination of the measurement the professional reads the data stored on the Flash Memory Card into his own PC for further analysis. The recorded data are presented as relative flow through the two nostrils.

Patient Groups

Patients > 12 years

The following patient groups can benefit from a rhinoflowmetry examination:

- Patients with temporary nasal complaints
- Patients with perceived obstructions in the nasal passage
- Patients with nasal dysfunctions
- Prospective OSA-patients
- Patients who are going through a functional rhinosurgery (before and after)
- Patients with a discrepancy between their clinical symptoms, endoscopic findings and diagnosis based on conventional nasal analysis.
- Patients suffering from allergy.